



Product datasheet

DCB-001NI

UV Disinfection Box



Features

Features of conventional products Safety instruction:

- 1. With UV tubes in two sides of the Box, 360 degree circulation disinfection
- 2. Can be a power supply for other devices by USB charging during disinfection process
- 3. Aromatherapy function to provide fragrance for the disinfection object
- 4. Indicator lights and voice prompts to ensure safe and correct use
- 5. Safe and portable with DC5V supply
- 6. Large capacity disinfection, applicable for up to 6 inch mobile phones, watches, toothbrushes, masks, underwear, jewelry and other small objects

Applications

Super sterilization and disinfection box



Good Choice As a Gift



Product datasheet

Technical data


Electrical data

Max capacity	6.5 inch mobile phone
Input voltage	DC 5V
Input current	2A
Power(diffuse)	1W
Power(sterilize)	2W
Max power	9W
Dimensions of unit	197*153*49MM
Dimensions of inner cavity	177.5*103.5*30MM
Ultraviolet wavelength	253.7nm

Directions for use

How to use:

Connect the power supply cable, put the object into the disinfection box, press the Sterilize button, then the object begins to be sterilized, and the Indicator lights is showing the disinfection time left. After 3 minutes, disinfection is completed and the UV tubes will turn off automatically. If you like aromatherapy, put essential oil in the aroma inlet, then press the Incense button anytime.




- Keep the unit far from flame, water and corrosive chemicals.
- Do not clean the device with organic solvent.
- Do not stare at the disinfection box for a long time.

Product datasheet

Dimensions & weight

Logistical Data

Package



Packaging unit (Pieces/Unit)	Dimensions	Gross weight	Volume
Small box 1	16*6*20.3cm	0.45Kg	0.002m ³
Outer Carton 30	49.5*42*33cm	14.5Kg	0.068m ³